


The Pioneering Medical Informatics  
Company of India...




[ [www.medsynaptic.com](http://www.medsynaptic.com) ]

**MED** < **SYNAPTIC**

# ScanDoc™ 3D

The easiest, fastest and best 3D software for Non DICOM / DICOM CT and MRI

- Compatible with any Non DICOM (Analog) or DICOM CT and MRI machine of any manufacturer
- Extremely fast 3D rendering
- Very easy to use
- Build 3D in 2 clicks
- MIP and MPR functions available
- Unique 3D slicer feature which allows user to move in the image volume and pinpoint the exact position of lesion
- Special algorithm to reduce opacity levels increasing quality of images
- Advanced smoothing algorithm to smooth the surface of the 3D images
- Apply color mapping to 3D images
- Rotate 3D images in all 3 axes X, Y, Z
- W/L, Zoom and Pan the 3D images
- Measurement function available
- Crop function available to remove unwanted data
- Print 3D images in various formats on paper
- Able to export the animated 3D volume to AVI format
- Able to save 3D images in jpeg format
- No special computer hardware required



47/6, "Vasantika", Opp Yena Bungalow,  
Paud Road, Kothrud,  
Pune - 411038, MS (India)

**Medsynaptic  
Private Limited**

Tel : 91 - 20 - 25443349 / 9881495351 / 52  
Fax : 91 - 20 - 25443913  
Email : sales@medsynaptic.com

[ w w w . m e d s y n a p t i c . c o m ]